




Visual Guide to iPhone Pairing




6. Ensure Bluetooth is turned on




5. Press "accessibility"




4. Press "hearing devices"



2. Wait for devices to populate in the search bar
 - a. Ensure aids have recently been restarted so they are able to be 'found'
 - b. If an old pairing is present at this step, click through and "forget" the old devices to continue pairing process



3. Tap the hearing aids under the "MFi Hearing Devices section once both are listed (R+L)



1. Press "pair" button when prompted (once for each device)
 - a. To then reconnect in any manufacturer app, just open the app and follow instructions provided at startup or check troubleshooting section for app-specific assistance